



# Non-invasive maternal blood test for identifying chromosomal causes of miscarriage



*Analysis of circulating cell-free fetal DNA by Next  
Generation Sequencing*



## POCADVANCE

cfDNA Karyo

[www.pocadvance.it](http://www.pocadvance.it)




## CLINICAL BACKGROUND: MISCARRIAGE AND DIAGNOSTIC RATIONALE

Miscarriage is the most frequent complication of the first trimester of pregnancy, with an estimated incidence of **15% to 25%** of clinically recognized pregnancies.

Although the etiological factors are multiple, approximately **60% of cases** are attributable to **chromosomal abnormalities** of the product of conception (POC).<sup>1–2</sup>

Identifying a chromosomal abnormality as the cause of pregnancy loss has significant clinical value: it enables **recurrence risk estimation** and may reveal **familial chromosomal rearrangements** that can predispose the couple to recurrent miscarriage or to the birth of offspring affected by congenital malformations and/or intellectual disability.

In this context, since 2016 the **American College of Obstetricians and Gynecologists (ACOG)** and the **Society for Maternal-Fetal Medicine** have recommended chromosomal analysis on placental tissue, amniotic fluid, or product of conception in all cases of intrauterine fetal demise or perinatal mortality.<sup>3</sup>



1. D'Ippolito S. et al., 2017

2. Rosenfeld J.A. et al., 2015

3. Committee Opinion No. 581, 2016



# POCADVANCE cfdNA Karyo

## THE POCADVANCE cfdNA KARYO TEST

**POCADVANCE cfdNA Karyo** is a non-invasive genetic test that, through the analysis of **circulating cell-free fetal DNA (cfdNA)** obtained from a maternal peripheral blood sample, allows determination of the **molecular karyotype** of the product of conception. The test uses high-sensitivity, high-resolution **massively parallel sequencing technologies (Next Generation Sequencing – NGS)** to detect fetal chromosomal abnormalities that may be responsible for the miscarriage event.

### Advantages of cfdNA analysis

- ✓ **Non-invasive:** performed with a simple maternal blood draw
- ✓ **Clinically useful timing:** can be performed before any uterine evacuation procedure
- ✓ **Applicable early:** potentially feasible in the **early stages** of pregnancy ( $\geq 5$  weeks)




## CIRCULATING CELL-FREE FETAL DNA (cfdNA): BIOLOGICAL RATIONALE

During pregnancy, the placenta physiologically releases DNA fragments into the maternal circulation (apoptosis), starting from approximately the **5th week of gestation**. This material is defined as circulating **cell-free fetal DNA (cfdNA)**.

Close to a miscarriage event, placental tissue may continue to release cfdNA into the maternal blood, making it possible to perform **non-invasive** investigations immediately after the ultrasound diagnosis of spontaneous pregnancy loss, with the aim of obtaining information useful on the potential genetic contribution.



## TECHNOLOGY: HIGH RESOLUTION AND ADVANCED BIOINFORMATICS





**POCADVANCE** cfDNA Karyo

## RESULT INTERPRETATION



### **POSITIVE**

**Chromosomal abnormality detected.**

A numerical and/or structural chromosomal abnormality has been identified in the fetal karyotype.



### **NEGATIVE**

**No chromosomal abnormality detected.**

No chromosomal abnormalities were identified; the fetal profile is compatible with a euploid karyotype.





# POCADVANCE

cfDNA Karyo

## INDICATIONS FOR TESTING

**POCADVANCE cfDNA Karyo** is particularly indicated in the following clinical settings:

- ✓ **spontaneous miscarriage** (anembryonic gestational sac, absence of cardiac activity, miscarriage in progress) when a timely, non-invasive investigation is desired;
- ✓ history of **recurrent miscarriage**
- ✓ **repeated pregnancy losses** in the absence of an exhaustive clinical explanation.



## HOW TO ORDER THE TEST



Kit  
request



Documentation  
completion



Sample  
collection



Sample  
shipment



Report  
delivery

### Blood draw timing

The blood sample should be collected **as soon as possible** after ultrasound confirmation of pregnancy loss and **before** any uterine evacuation procedure; alternatively, **within 24 hours** of expulsion of the products of conception.

The test can be performed in **singleton or twin pregnancies** (monozygotic or dizygotic) at **≥ 5 weeks** of gestation.

Turnaround time



2-3 working days

GENOMICA is a highly innovative company with extensive technical and scientific expertise, active in both clinical applications and research. Supported by a team with over 20 years of experience in molecular diagnostics, GENOMICA combines cutting-edge technology with a strong commitment to innovation, delivering increasingly accurate and accessible diagnostic services.



Over **100.000 genetic tests/year**



Laboratories with **groundbreaking technologies** and high quality standards



Dedicated **R&D Team**



International **Partnerships**



**Personalized genetic counseling** with genetic counselors experts in discussing genetic test results and familial risks



**20+ years experience** in prenatal molecular diagnostics

## LABORATORIES

Rome: Via Arduino 38 - 00162 Tel.: 06.21115020  
E-mail: [info@genomicalab.it](mailto:info@genomicalab.it)  
[www.genomicalab.it](http://www.genomicalab.it)

## REGISTERED OFFICE

Rome: Via Arduino 38 - 00162  
PEC: [info@pec.genomicalab.it](mailto:info@pec.genomicalab.it)  
VAT No.: 14554101007 - REA: RM - 1530210

Visit the website  
dedicated to the test

